



Disposables Hand Switch Pencils (Autoclave)

Electrosurgical Pencil:	Push Button	Single-Use:	Disposable
Hand Switch:	Handle Stainless Steel	Electrode:	3m Cable

Description:-

- 1) Model: HT-2A
- 2) Length of handle: 170mm
- 3) Length of blade: 4mm
- 4) Diameter of blade: 2.3mm
- 5) Length of wire: 2.8m

Specifications:

- 1) Max. output current: 1A
- 2) Frequency confine: 0-1.0MHz
- 3) Current resistance: 200
- 4) Ways of aseptic technique: Sterile or irradiation
- 5) Term of sterility: 2 years

Packing:

100pcs/carton, Carton Dimensions: 46 x 26 x 28cm, N.W.: 7.7kg, G.W.: 8.7kg. Standard Push Button Electrocautery Pencil with Electrode, including Needle, Blade, Ball, Loop or PTFE Non-Stick Electrode, 3m Cable with one Pin or 3 Pin Plug.

It is compatible with most of generators like Valleylab, Erbe, Bowa, Martin, Bovie, etc.
Standard Push Button Electrocautery Pencil with Electrode, compatible for kinds of units.

CE MARK + ISO-9001-2000 CERTIFIED

Address

Mustafa Town, Fateh Garh, S.I.E Sialkot - 51310, Pakistan

E-Mail: info@g14surgical.com | g14surgical@brain.net.pk | g14surgical@gmail.com Web: <http://g14surgical.com>

Tel: +92- (0) - 52- 3 256108, 3258365 Fax: +92-(0)-52-3 240858, Cell: +92-(0)-345-6724741