



G14 SURGICAL SUPPLY

Manufacturers & Exporters of all kinds of Surgical, Dental and Laparoscopy Instruments

LAPAROSCOPIC CLIP APPLIERS



Features for 5 & 10 MM Ø CLIP APPLIERS:-

1. Modular, 3-piece screw-apart instrument for simple cleaning & maintenance.
2. Constructed of the highest grade stainless steel.
3. Durable, non-deforming jaw for consistent clip application.
4. Complete assembly for easy cleaning and sterilization.
5. Precisely manufactured to insure highest quality standards and longevity.
6. Full hand grips for safe and improved instrument handling.
7. Available in three sizes - small, medium, and medium large.



1. Clip-Applicator for Medium clips, Rotatable.
2. Clip-Applicator for Medium large Clips, Not Rotatable.
3. Clip-Applicator for LT 300 clips, Not Rotatable.

Product Code	Product Code	Clip Size/No. of Clips	Clip Dimensions		Overall Clip Applier Length
			Open	Closed	
	G14 - # 1099	Small 20	2.1 mm	3.8 mm	23.8 cm
	G14 - # 1100	Medium 20	4.3 mm	6.0 mm	23.8 cm
	G14 - # 1101	Medium 20	4.3 mm	6.0 mm	29.2 cm
	G14 - # 1102	Medium 30	4.3 mm	6.0 mm	29.2 cm
	G14 - # 1103	Large 20	6.3 mm	10.8 mm	33.7 cm

CE MARK + ISO-9001-2000 CERTIFIED

Address

Mustafa Town, Fateh Garh, S.I.E Sialkot - 51310, Pakistan

E-Mail: info@g14surgical.com | g14surgical@brain.net.pk | g14surgical@gmail.com Web: <http://g14surgical.com>

Tel: +92- (0) - 52- 3 256108, 3258365 Fax: +92-(0)-52-3 240858, Cell: +92-(0)-345-6724741