



G14 SURGICAL SUPPLY

Manufacturers & Exporters of all kinds of Surgical, Dental and Laparoscopy Instruments


LAPAROSCOPIC ELECTRODES

Features Specifications: G14 Laparoscopic Electrodes

Tungsten loop: Stainless Blade

-Electrodes: Disposable

-Electrodes including tungsten loop, triangle electrode, square loop electrode, needle, ball and standard blade. 2.36mm standard diameter shaft, 5mm is optional also.



Electrodes including Loop, Blade, Needle, Ball, etc, compatible with all standard Electrosurgical Pencils.

CE MARK + ISO-9001-2000 CERTIFIED

Address

Mustafa Town, Fateh Garh, S.I.E Sialkot - 51310, Pakistan

E-Mail: info@g14surgical.com | g14surgical@brain.net.pk | g14surgical@gmail.com Web: <http://g14surgical.com>

Tel: +92- (0) - 52- 3 256108, 3258365 Fax: +92-(0)-52-3 240858, Cell: +92-(0)-345-6724741