





LAPAROSCOPIC TC-NEEDLE HOLDER



TC-Needle Holders

With TC Super grip inserts
 Axial handle with closing and reset spring
 Online ring handle with synchronized mechanism
 Designed with TC-inserts for long life
 Ergonomic and user-friendly
 Available in three different jaw designs
 Luer Lock irrigation adapter

5mm Precision Laparoscopic Needle Holders, with Tungsten Carbide inserts, Made in Pakistan, 33 and 45cm lengths.

	Needle Holder, Solid version with Axial Handles, Short 8 mm	G14 # 1105
	Needle Holder, Cleaning Port, Break proof, Axial Handles, Long 10 mm	G14 # 1106
	Needle Holder, Cleaning Port, Break proof, Axial Handles	G14 # 1107
	Needle Holder, Cleaning Port, Break Proof, Axial Handles	G14 # 1108

CE MARK + ISO-9001-2000 CERTIFIED

Address

Mustafa Town, Fateh Garh, S.I.E Sialkot - 51310, Pakistan

E-Mail: info@g14surgical.com | g14surgical@brain.net.pk | g14surgical@gmail.com Web: <http://g14surgical.com>

Tel: +92- (0) - 52- 3 256108, 3258365 Fax: +92-(0)-52-3 240858, Cell: +92-(0)-345-6724741