



G14 SURGICAL SUPPLY

Manufacturers & Exporters of all kinds of Surgical, Dental and Laparoscopy Instruments

LAPAROSCOPIC TROCARS

G14 Surgical Laparoscopic Trocars are in general made from Stainless Steel & metal, Reusable and disposable Trocars are constructed by a combination of Stainless Steel & metal. The tip of disposable Trocar has a two edged blade. these are very effective at penetrating the abdominal wall by cutting the tissue as they pass through. Most of the disposable metal trocar has a spring loaded mechanism that withdraws the sharp tip immediately after it passes through the abdominal wall to reduce the incidence of injury of viscera.

Trocars of different sizes and diameter depending upon the instrument for which it is used. The diameter of Trocar range from 5-MM to 10-MM. The metal Trocar has different type of tips i.e. Pyramidal Tip, Eccentric Tip, Conical Tip or Blunt Tip depending on the surgeon's experience.



Size: 10mm

Made of Stainless Steel & Metal Trocars conform to the specifications of Laparoscopy surgery. All these laparoscopy trocar provide smooth and consistent cutting pressure with the cutting edges and piercing points.

All our laparoscopy trocar are designed for single-use & re-useable only. Design & specifications may vary as per the clients' requirements. The dimensional and material specifications are given below:



Size: 5mm

Address

Mustafa Town, Fateh Garh, S.I.E Sialkot - 51310, Pakistan

E-Mail: info@g14surgical.com | g14surgical@brain.net.pk | g14surgical@gmail.com Web: <http://g14surgical.com>

Tel: +92- (0) - 52- 3 256108, 3258365 Fax: +92-(0)-52-3 240858, Cell: +92-(0)-345-6724741